



PLANMECA
ProMax 3D s



Profound understanding of anatomy

The unique Planmeca ProMax 3D product family offers equipment for all maxillofacial imaging. All volume sizes from the smallest special cases to whole head images are available. Planmeca ProMax 3D s is designed to obtain complete information on patient anatomy in the minutest detail. The unit provides 3D, digital panoramic and cephalometric imaging, as well as advanced imaging software tools to comply with every possible need in dental radiology.





The amorphous silicon flat panel produces accurate, distortion-free images for 3D reconstruction. Unlike image intensifier sensors that use old vacuum tube technology and multi-step focusing, flat panels use single step image readouts with no geometric distortion, no loss of sensitivity, and therefore no need for frequent calibration.

Planmeca's proprietary 3D reconstruction algorithm converts the original 2D transillumination images to 3D volume study and is the core component for high quality

3D imaging. It handles high contrast objects, like amalgam fillings, in a special way in order to produce undisturbed study views.

The reconstructed image volume consists of more than 200 million voxels. These voxels are isotropic, which enables accurate 1:1 measurements and ensures geometric relations throughout the image. The extremely small voxel size, $100^3 \mu\text{m}^3$, provides a detailed high-resolution 5 lp/mm image without artefacts.



In modern dentistry, the demand for surgical implant treatments is steadily growing. The growth trend has created a need for a more advanced X-ray imaging device. Planmeca ProMax 3D s, a Cone Beam Volumetric Tomography (CBVT) unit, is expressly designed to comply with the needs of modern surgical dentistry. It supplies clear, dependable imaging in a three-dimensional format with limited patient radiation dose.

Thanks to its small footprint, Planmeca ProMax 3D s makes effective three-dimensional imaging possible in every dental practice. Planmeca ProMax 3D s is a genuine all-in-one unit including digital panoramic, digital cephalometric, and 3D digital imaging, all in the same unit, saving office space and investment costs. This innovative, versatile, and dynamic imaging

device will open up new possibilities for dentists on-site in their offices.


Planmeca ProMax 3D s uses new Cone Beam Volumetric Tomography technology, in which the X-ray beam used is shaped in the form of a cone or a pyramid. The CBVT technology takes the whole volume needed in a single semicircle scan, as opposed to a medical CT that takes multiple axial slices in multiple full circle scans.

During the scan, each image is made using a short X-ray pulse instead of continuous radiation. The total scanning time is 18 seconds for one volume, but the actual exposure time is only 3 seconds at shortest. This technique reduces patient radiation dose considerably and forms stroboscopic X-ray effect which, together with the shortened rotation scan


(only 200 degrees), virtually eliminates artefacts contributing to outstanding image quality.

The Planmeca ProMax platform's unique SCARA technology (Selectively Compliant Articulated Robot Arm) enables free formation of the imaging geometry. Planmeca's patented, computer-controlled SCARA robotic arm can produce any movement pattern required, ensuring perfectly accurate and reliable image volume positioning and sizing. All controls are made on a full colour graphical user interface in the language of your choice.


Thanks to the original, technologically advanced design, any Planmeca ProMax can be upgraded into a 3D Cone Beam Volumetric Tomography unit.




Wisdom tooth extraction
The image reveals that the extraction would be difficult, as the mandibular canal is located lingually to the roots.




TMJ study
The condyle is displayed sharply. The image clearly shows the condition of the temporomandibular joint. Malign finding can be seen in the right condyle head.




Impacted canine
An impacted maxillary left canine found behind the incisors.



Supernumerary tooth
An extra ninth molar found in maxilla.



Implant case
The upper right first molar is missing. The image clearly shows that there is just enough bone to place an implant.



Impacted canine
An impacted maxillary left canine found behind the incisors.

Unequaled imaging programs

Planmeca ProMax 3D s complies with a multitude of diagnostic requirements: those of implantology, endodontics, periodontics, orthodontics, dental and maxillofacial surgery, and TMJ analysis.

The 50 x 80 mm image size is optimum for most implantology applications. The basic volumes can also be stitched together to generate an image up to 75 mm in width.

Planmeca ProMax 3D s produces high-resolution volumetric studies of the mandible and maxilla for analysing the bone structure available, the location of the mandibular canal, and the correct position

for the implant. Pre-surgical planning will reach a new level of precision, as the prospective site is visible in all three imaging planes: sagittal, axial, and coronal.


Third molars, maxillary cuspids, supernumerary teeth, and impactions challenge the clinician to identify the tooth's orientation. With Planmeca ProMax 3D s all anatomical structures, angles and orientations become clearly visible.

Planmeca ProMax 3D s studies accompanied by digital cephalometric images provide full visualisation of all classes of orthodontic malocclusion. This is highly advantageous for orthodontic planning,


as time is saved and patient radiation dose reduced. Planmeca ProMax 3D s provides the image data in the correct anatomic 1:1 ratio, without need to correct for geometric magnification.

Planmeca ProMax 3D s also provides high-resolution TMJ studies for true and accurate evaluations of the joint arthritides, condylar morphology, and the condyle-fossa relationship.


With its high resolution (5 lp/mm) and advanced reconstruction technology, Planmeca ProMax 3D s establishes the new standard for 3D dental radiology.




Finding in canine area
A suspicious mass found behind the lower right canine.



Finding in sinus
Swelling in sinus mucosa.



Lesion in right mandible
A radiolucent lesion is clearly visible in the right mandible.
The probable cause is an unsuccessful root treatment.



Planmeca Romexis 3D Implant module

Planmeca Romexis 3D Explorer, the 3D image acquisition software for Planmeca ProMax 3D s, enables flexible multiplanar viewing in all three relevant projections: axial, coronal, and sagittal. The software incorporates a re-slicing feature, which enhances the projections and enables real-time three-dimensional viewing in the desired angle. A rendered 3D view provides a realistic overview of the anatomy.

With the Planmeca Romexis 3D Explorer software, each patient study can be stored on a CD with Planmeca Romexis 3D Viewer for others to view.

The Planmeca Romexis 3D Cross Sections module, which is available as an option, produces cross-sectional

images of anatomy along with the defined panoramic curve. The number of the cross-sections, axial views, and lateral views, and their exact positions can be freely chosen. The 3D Cross Sections module also includes the reconstructed panoramic view. The panoramic image is created from the acquired volume of data without the undesired artefacts, common in normal panoramic images. As the image is reconstructed using software, the user can determine the location and thickness of the focal trough.

The optional Planmeca Romexis 3D Implant Planning module offers tools for implant placing and nerve drawing. The implant placements are determined

with the help of realistic implant models of the most common implant brands. A drawing tool allows clear marking of the mandibular nerve.

Planmeca Romexis software has optional DICOM functionality, which allows 3D studies to be transferred to other implant planning software or any other software that receives images in DICOM format. Studies can also be transferred to PACS or to a high quality DICOM printer in the network.

Planmeca Romexis is a pure Java based software that runs in various operating systems and modern web environments.

Planmeca Romexis software

Planmeca Romexis is a complete dental imaging software, including all dental imaging modalities: intraoral, panoramic, cephalometric, 3D imaging, dental tomography as well as intraoral video and still camera imaging. With a complete set of tools for image viewing, enhancement, measurements, and annotations, Planmeca Romexis also improves the diagnostic value of radiographs. Printing, image import and export, and DICOM functionalities are also included.

Planmeca Romexis platform fully integrates digital imaging with the patient's other clinical data. The system provides direct image capture from Planmeca's X-ray equipment, and interfaces with 3rd party devices via TWAIN. Together with Planmeca's X-ray equipment, Planmeca Romexis provides a unique safety feature especially useful for teaching environment: the X-ray image capture is inhibited until the supervisor has approved the student's image capture request.

Planmeca Romexis computer recommendations

	Planmeca Romexis client work station	Planmeca Romexis server
Processor	2 GHz Core Duo or equivalent	3 GHz Core Duo or equivalent
RAM	3 GB	3 GB
Hard disk space	40 GB	2 x 500 GB (RAID1 mirroring)
Graphics card	ATI or NVIDIA, 128 MB minimum memory	Not required
Monitor	1280 x 1024	1024 x 768
Peripherals	CD R/W or DVD R/W drive	CD R/W or DVD R/W drive
Backup medium	None necessary	DAT or equivalent
Operating system	Windows XP, Windows 2003, Windows Vista, Mac OS X, Linux Mac OS X / Linux support subject to contract	Windows XP Pro, Windows 2003 Server, Windows Vista
Other	Java platform (Java Virtual Machine 1.6 or later)	Java platform (Java Virtual Machine 1.6 or later)

The disk space requirements are dominated by digital images. Thus the space requirements vary, but a rough estimate is in the order of 1 MB per 2D X-ray image, 7–9 MB per extraoral image, depending on a variety of image specific factors, and 250 MB per 3D image.

It is recommended to use the same computer for application server and for database server. If Planmeca Romexis server computer is also used for client activities, the hardware should meet both client and server specifications.

Values presented above are recommended minimum specifications, not meeting them may in some cases cause degraded performance.


DICOM compatibility

- Media Storage – saving images into removable DICOM media
- Print – printing images on film or paper with a DICOM medical printer
- Storage – saving images into DICOM image archive
- Query/ Retrieve – importing digital images from DICOM image archive
- Worklist – importing a patient list from DICOM patient management
- Storage Commitment – confirmation of a successful image storage


Planmeca ProMax 3D s

X-ray beam	Cone
Focal spot	0.5 mm, fixed anode
Image detector	Amorphous silicon flat panel
Gray scale	15 bit
Detector resolution	630 x 1024 pixels, pixel size 127 µm x 127 µm
Voxel size	100 x 100 x 100 µm, isotropic 200 x 200 x 200 µm, isotropic
Image acquisition	Single 200 degree rotation
Total scan time	18 s, pulsed X-ray
Reconstruction time	30–150 s
Standard volumes (diam. x height)	Ø50 x 80 mm (child mode Ø42 x 68 mm) Ø50 x 50 mm (child mode Ø42 x 42 mm)
Stitched volume (diam. x height)	Ø100 x 110 mm (child mode Ø85 x 90 mm)
3D reconstruction server	Proprietary Feldkamp type back projection reconstruction algorithm Improved Artifact Removal (IAR) High Contrast Object Compensation (HCOC)

Example installation



Dimensions and space requirements



		Planmeca ProMax 3D s	Planmeca ProMax 3D s with cephalostat
Physical space requirements	Width	96 cm (38 in.)	194 cm (76 in.)
	Depth	125 cm (49 in.)	125 cm (49 in.)
	Height*	153–243 cm (60–96 in.)	153–243 cm (60–96 in.)
Minimum operational space requirements	Width	150 cm (59 in.)	215 cm (85 in.)
	Depth	163 cm (64 in.)	163 cm (64 in.)
	Height*	243 cm (96 in.)	243 cm (96 in.)
Weight		113 kg (lbs 248)	128 kg (lbs 282)

*The maximum height of the unit can be adjusted for offices with limited ceiling space.



Planmeca Oy designs and manufactures a full line of high technology dental equipment, including dental care units, panoramic and intraoral X-ray units, and digital imaging products. Planmeca Oy, the parent company of the Finnish Planmeca Group, is strongly committed to R&D, and is the largest privately held company in the field.



Planmeca Oy
 Asentajankatu 6 | 00880 Helsinki | Finland
 tel. +358 20 7795 500 | fax +358 20 7795 555
 sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations. Some products displayed above may not be available in all countries or areas. Rights for changes reserved.